TRAINING & WELLNESS Equine Nano Treatment System

Nano technology - treatment of the future!


The Equine Nano Treatment System is a revolution in Equine health! Treatment with minimal doses and maximimal results is possible with this system. This system surpasses all traditional treatments and focuses on the future: nano-technology.

Nanoparticles are so small that they are not visible to the human eye. One particle is thirty times smaller than a human skin cell. Nanoparticles are incorporated seamlessly by the skin because of their ability to pass through the horse's skin barrier. The particles are also taken up easily into the horse's respiratory system. Tablets, powders, syringes and ointments are no longer required for a successful treatment.

Skin and respiratory problems are a common phenomenon in the equine species. Often, these problems can only be partially cured. Orally administered drugs are largely lost in the digestive system, which is why big amounts of medicine need to be administered which can cause unpleasant side-effects. The Equine Nano Treatment System avoids the very senitive digestive system of horses which minimizes side-effects.



COLUMBUS



Wellness

The Equine Nano Treatment System can also be used for improving the horse's health. There are ampoules available with the following substances:

- Vitamins A, E and C
- Dead sea mineral molecules
- Hyaluronic
- Acid (HA)

Healing

In the near future all traditional treatments are supposed to my available for nano-therapy. On request we can provide you with information about the substances and active agents which are already available.

Columbus Professional Horse Equipment Eigentümer: Andreas Wegmann Lissendorfer Strasse 19 D - 54584 Gönnersdorf

All changes, errors and misprints reserved

Tel: 0049 (0)6597 961092 Fax: 0049 (0)6597 961093 Email: info@columbus-de.com Web: www.columbus-de.com

TRAINING & WELLNESS Equine Nano Treatment System

COLUMBUS



Use of the Equine Nano Treatment System

The Equine Nano Treatment System is very easy to use. The system consists of an airtight cabin which can be customized and built-in on request. The system is attached to the cabin with an airtight tube. The ampoules with the healing substances are inserted into the system. The Equine Nano Treatment System is delivered with a comprehensive guide. Applications, treatments, treatment times and periods are advised.

In addition to the airtight cabin we can also deliver a special inhalator for resperatory ailments which can be attached to the ENTS.

Why Equine Nano Treatment System?

The Equine Nano Treatment System is an entirely new product which provides new ways of wellness and medicine administration. This innovative technique offers new possibilities to ensure health and cure diseases. A scientific study regarding the effects of treating horses with bronchitis using the Equine Nano Treatment System has demonstrated the following benefits:

1. Minimal dosage, maximal results

Only a small amount (on average 1/5) of the traditionally used dosage of medicine is necessary to get the same effect. This means faster and more efficient results as well as minimized side-effects.

2. Pain free treatment

Treatment using the Equine Nano Treatment System is pain free for your horse. Traditionally medicines are administered intravenous (directly into the vein) or intramuscular (directly into the muscle), which is often painful for the horse and can lead to dangerous situations. Treatment using the Equine Nano Treatment System is therefore not only pain free and safe, but also minimizes stress.

3. Easy treatment

Using the Equine Nano Treatment System, You as owner can treat your horse yourself. Traditional administering often need to be performed by a vet. The use of the Equine Nano Treatment System is simple: Your vet prescribes the medicine, you treat your horse.

4. Safety and doping

Medicines are supplied to you in ampoules. These ampoules are filled at the University of Graz. They are the only ones who perform a nano-aerosol filling on sterile basis. The dosage cannot be changes after the filling in Graz. Therefore, it is reliable and accurate.

Columbus Professional Horse Equipment Eigentümer: Andreas Wegmann Lissendorfer Strasse 19 D - 54584 Gönnersdorf

All changes, errors and misprints reserved

Tel: 0049 (0)6597 961092 Fax: 0049 (0)6597 961093 Email: info@columbus-de.com Web: www.columbus-de.com

TRAINING & WELLNESS Equine Nano Treatment System

COLJMBUS

ENTS – a revolution for equine health!

- \checkmark effective health improvement of your horse
- \checkmark treatments are individualized for each horse
- ✓ easy to operate with touch screen
- \checkmark built-in ventilation system of high quality
- ✓ impact resistant materials
- ✓ easy to install
- ✓ low on maintenance
- ✓ several safety measures
- ✓ all steel parts are galvanized according to NEN-EN-ISO 1461-99





Produktdetails ENTS


Technical Data


| Connection voltage | 1PH + N + PE 230 V 50 Hz |
|----------------------------|--------------------------|
| Circuit breaker: | 16 A |
| Power consumption | 0,25 kW |
| Noise level at distance of | |
| 4m 79/113/EWG | 31 dB |
| | |

Materials and dimensions

| Material | steel, galvanized according to |
|------------------------------|--------------------------------|
| | NEN-EN-ISO 1461-99 |
| Powder coating | RAL9010 - white |
| Length of cabin | 350cm |
| Width of cabin | 210cm |
| Height of cabin | 250cm |
| Amount of doors | 4pcs |
| Windows made of safety- | |
| glass | \checkmark |
| Rubber flooring | \checkmark |
| Weight | 1000kg |
| Equipment and Extres | |
| Equipment and Extras | |
| Cable and Plug | |
| Control box with touch scree | |
| Horse crush | |
| Ventilation system | \checkmark |
| CE certified | \checkmark |
| Powder coating in | |
| preferred colour | 0 |

Optional





Columbus Professional Horse Equipment Eigentümer: Andreas Wegmann Lissendorfer Strasse 19 D - 54584 Gönnersdorf

 $\sqrt{}$

Inklusive

All changes, errors and misprints reserved

0

Tel: 0049 (0)6597 961092 Fax: 0049 (0)6597 961093 Email: info@columbus-de.com Web: www.columbus-de.com